
5/98-Sf-

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This unit can amplify DC voltages and AC voltages of any curve form up to 1000 times to a maximum output voltage of 15 V ( $\approx 30 V_{pp}$ ), for a maximum load of 2 A. A built-in 5-V DC voltage source means that no additional voltage source is required for the input-side connection of e.g. the carbon microphone (586 28) or a logic circuit.

The AC/DC-amplifier is short circuit-proof and is equipped with an automatic overheating cut-out and separate LEDs to indicate when the system is overdriving and the maximum output current is exceeded, to ensure a high degree of protection against damage due to improper use. An additional LED indicator showing when the output voltage is zero permits offset compensation without connecting an external meter.

#### Possible applications:

- Amplification of slight voltage changes in the audio frequency range (e.g. microphone or oscillator circuit).
- Measuring amplifier for DC and AC voltage in the millivolt range.
- Power amplifier for function generators (e.g. for generating low-frequency alternating magnetic fields with currents up to 2 A).

#### 1 Safety notes

115 V AC : s. cat.no. 532 20 NA.

Compensate the offset before connecting a loudspeaker (cf. section 3.1).

Do not obstruct the flow of air to the cooling fins of the heat sink during operation.

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## 2 Description, technical data (Fig. 1)

- ① Power switch
  - ② Pair of 4-mm input sockets (INPUT), bottom socket connected to output ground, floating  
Input impedance: 1 M $\Omega$   
Input protected up to 230 V<sub>rms</sub>
  - ③ AC/DC-selection switch:  
Position DC: DC voltage amplification  
Position AC: Amplification of the AC component of the input voltage  
Frequency range: 1.6 Hz ... 100 kHz (-3 dB)
  - ④ Continuous attenuator 0 ... 100 % for the gain factor set with ⑤; right stop calibrated
  - ⑤ Gain selection switch  
Gain factors: x 1, x 10, x 100, x 1000, corresponding to input voltage ranges from 15 V ... 15 mV  
Accuracy: 1 % for attenuator at full right stop ④
  - ⑥ Offset potentiometer for compensating the offset voltage. LED indicator ⑨ (see also section 3.1)
  - ⑦ Pair of 4-mm output sockets (OUTPUT), bottom socket connected to input ground socket; floating  
Output voltage: maximum  $\pm 15$  V ( $\cong 30$  V<sub>pp</sub>)  
Output current: maximum 2 A ( $\cong 4$  A<sub>pp</sub>), electronically limited  
Output impedance: 50 m $\Omega$  + 5.6 mH  
Output short-circuit proof, proof against external voltages up to 15 V<sub>rms</sub>  
Noise 0 ... 20 kHz < 20  $\mu$ V<sub>pp</sub>, with reference to input
  - ⑧ Red LED:  
Lights up for output currents  $I_a > 2$  A (current limiting activated),  
flashes to indicate overheating cut-out; cf. section 3.2
  - ⑨ Green LED:  
Lights up for output voltages  $|U_a| < 50$  mV; cf. section 3.1
  - ⑩ Red LED:  
Lights up for overdriving of amplifier ( $|U_a| > 15$  V)
  - ⑪ 4-mm socket with bias voltage approx. 5 V DC via 18 k $\Omega$ ; cf. section 3.3
- Power consumption: 85 VA  
Mains voltages: 230 V AC/115 V AC  
Dimensions: 20 cm x 14 cm x 23 cm  
Weight: approx. 4 kg

Appliance plug with integrated holder for primary and spare fuse on rear of housing. See fuse protection rating plate above on rear of device.

The amplifier is supplied with fuses for mains voltage 230 V AC installed.

Amplifier supplied complete with: mains cable, 2 fuses for mains voltage 115 V AC.

## 3 Operation

### 3.1 Setting up the unit

Connect the input ② and output ⑦ as appropriate for the experiment; however, only connect a loudspeaker after carrying out offset compensation.

Set the AC/DC-selector switch ③ to DC amplification ("DC") or AC amplification ("AC").

Use the gain selector switch ⑤ to set the gain so that LED ⑩ does not light up during operation. If necessary, reduce the gain using the attenuator ④.

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When using the device as a measuring amplifier, set the attenuator ④ to the right stop (position "Cal"), so that the setup has exactly the gain factor set with ⑤.

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Offset compensation:

Every amplifier has an output voltage  $U_a$  which is not equal to zero, even when no input signal is present. Its DC voltage component is referred to as the "offset"; the AC voltage component usually consists of a noise signal, sometimes with a hum voltage superposed on it. Using the potentiometer ⑥, you can compensate the offset, i.e. reduce the output voltage to approximately zero.

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Short-circuit input ② or set attenuator ④ to zero (left stop). Adjust the offset potentiometer ⑥ so that the green LED ⑨ lights up (indicating an output voltage  $|U_a| < 50$  mV). Readjust the offset compensation each time after switching the gain factor (switch ⑤).

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### 3.2 Warning indications of red LEDs ⑧ and ⑩

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Display	Cause	Remedy
LED ⑩ lights up	Output overdriven: $ U_a  > 15$ V	Turn attenuator ④ or gain selection switch ⑤ to left
LED ⑧ lights up	Output current $> 2$ A; current limiting possibly actuated	Turn attenuator ④ or gain selection switch ⑤ to left or increase load resistance
LED ⑧ flashes	Internal device temperature too high; current automatically limited to approx. 5 mA	Improve cooling; increase load resistance

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When the device is cool, current limiting cuts in at approx. 2.5 A. The amplifier can then be loaded up to this level with no distortion, even when LED ⑧ lights up.

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Illumination of ⑧ when output ⑦ is unconnected is irrelevant

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#### 3.3.1 Connecting a carbon microphone (Fig. 2)

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Connect the microphone to the pair of input sockets ②, and connect the top input socket with the 5-V socket ⑪. Set the AC/DC selector switch ③ to "AC" and the gain to x 100.

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Connect digital circuits with open-collector output to input sockets ② and connect the 5-V socket ⑪ to the top input socket. Set the AC/DC selector switch ③ to "DC" and the gain to x 1.

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#### 3.3.3 Using the amplifier as a power supply 15 V/2 A

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Connect the 5-V socket ⑪ to the top socket of input ② and set AC/DC-selector switch ③ to DC; set the gain to x 1, for output voltages above 5 V (max. 15 V) set the gain to x 10. Roughly adjust the output voltage with the attenuator ④, then carry out the fine adjustment with the offset potentiometer ⑥.

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### 3.4 Experiment examples

See Fig. 2 to Fig. 4 on pp. 5/6

## 4 Replacing the fuse and converting the mains voltage

### 4.1 Replacing the primary fuse

- Pry out insert **a** with holder for primary fuse **b** and spare fuse **c**) (Fig. 5.1).
- Replace the defective fuse **b** using a new fuse **c** which has been checked for the correct rating (Fig. 5.2).
- Insert the new spare fuse **c** and put the insert **a** back in the device.

### 4.2 Converting the amplifier for mains voltage 115 V AC

*Caution! Unplug the amplifier from the mains!*

Unscrew the housing screws **e** on the bottom of the device using a Phillips head screwdriver (size 2) (Fig. 6.1).

Place the device on the bench right side up and remove the top section of the housing **f** (Fig. 6.2).

Take the PC board in the 230 V **h** configuration out of the blue plug-in socket **g** at the transformer (Fig. 6.3).

Turn the board **h** so that the printed value 110 V (for 115 V mains voltage) is at the bottom left above the 1 on the side of the plug-in socket **g** (Fig. 6.4).

Push in the board **h** and reclose the housing.

Do not forget to insert the appropriate primary fuse for the new mains voltage (see technical data).

See section 4.1 (Fig. 6.1 and 6.2) for fuse replacement instructions.


Fig. 5.1


Fig. 5.2



Fig. 6. 1 - 6.4


**Fig. 2 - Fig. 4 Versuchsbeispiele zu 522 61**  
**Examples of Experiments for 522 61**  
**Exemples d'expériences pour 522 61**  
**Ejemplos de experimentos por 522 61**

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Model experiment for piezoelectric microphone: generation of an AC voltage from the change in the pressure on a Rochelle salt crystal (587 25) produced by an oscillating tuning fork. Amplification via AC/DC-amplifier for sound reproduction using a loudspeaker (587 08).

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Connection of a carbon microphone (586 28) to AC/DC amplifier: 5 V socket (a) for microphone voltage supply, connected to signal input (b).

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Application as a power amplifier for acoustically demonstrating the Barkhausen effect; magnetization and demagnetization of soft iron (binding wire ③) by a slowly changing magnetic field in coil ⑥; power amplification of the voltage surges induced by changes in the magnetization state of the soft iron in coil ④ due to the Barkhausen jumps. The rising and falling loudspeaker noise is an acoustic indication of the Barkhausen effect  
 Utilisation comme amplificateur de puissance pour démontrer l'effet Barkhausen: magnétisation et démagnétisation du fer doux (fil pour fleurs ③) par la lente variation du champ magnétique de la bobine ⑥; amplification de puissance des chocs de tension induits dans la bobine ④ à la suite des sauts de Barkhausen au changement de la magnétisation du fer doux; l'augmentation et la diminution du bruit dans le haut-parleur caractérisent les sauts de Barkhausen.

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