Fusenet - WP7.A Documentation  
Hands-on Experiment: Magnetic Diagnostic

Edition 1.2

Eindhoven University of Technology <sup>1</sup>

23 April 2013

<sup>1</sup>Contact Information: [j.w.oosterbeek@tue.nl](mailto:j.w.oosterbeek@tue.nl)

# Contents

# 1 Introduction

To support the lectures in the Master track "Science and Technology of Nuclear Fusion" at the Eindhoven University of Technology, a PlasmaLab is in operation. In the lab the students experiment with the concepts that are introduced in the courses and are able to develop their experimental skills. The laboratory is a joint effort of all plasma groups in the department of Applied Physics. A total of 10 experiments are available:

1. Waves in Plasma
2. Zeeman Splitting
3. **Magnetic Diagnostic**
4. Spectroscopy
5. Paschen Curve
6. Probes
7. Interferometer
8. Laser Spectroscopy
9. Resonance Spectroscopy
10. Fusor

The first three of these experiments are co-funded by the the Fusenet Association. This report provides technical documentation on the **Magnetic Diagnostic** experiment such that Fusenet Association partners can build a duplicate.

## 1.1 Topic and teaching aims

Magnetic field measurements provide important machine parameters in tokamaks such as the plasma current and the location of the plasma. Local fluctuations in current density, such as MHD instabilities and sawteeth, are extensively probed with magnetic pick-up coils. Magnetic diagnostics have been around for a long time but new challenges also appear. For instance, an accurate measurement of the plasma current becomes more difficult as the pulse length increases due to drift and interference that add up in the integrators. Another difficulty is where to best mount the pick-up coils: inside the vacuum vessel is preferred to avoid losses due to Eddy currents and attenuation in the vessel wall, but inside the vessel the environment is very harsh and access in case of repairs is limited.

In this experiment the fundamentals of magnetic diagnostics can be analysed. In the experiment a plasma current is mimicked using a current carrying conductor. The current in the wire can be gradually ramped up, kept constant, and ramped down again resembling a plasma pulse. The wire runs perpendicular to a plate with an array of magnetic pick-up coils. The wire may also be moved along a slit in the plate. The signals of the probes can be displayed directly on an oscilloscope, can be amplified and integrated with discrete electronics, or can be digitised, and processed with software. In the experiment the students are asked to i) calculate the current through the wire given the number of windings in the pick-up coils, ii) give an expression for the location of the wire using two coils on opposite sides of a slit in the plate through which the wire is moved, and iii) give an expression for the location of the wire using all available pick-up coils. Item iii) may be evaluated experimentally too in case the students have sufficient experience programming in LabView<sup>®</sup> to program this expression. A second exercise, but not yet fully implemented at the time of writing the report, is a measurement with a very fast probe ( $\mu\text{s}$ -scale) to assess Eddy currents and attenuation in a stainless steel tube that can optionally be mounted around the wire. A very fast change in the current in the wire is simply generated by switching the power supply on and off. Finally the students are asked to comment on design considerations for the slow measurement system (plasma current) and the fast measurement system (MHD-fluctuations).

## 1.2 Overview of experiment

The overall experiment is depicted in Fig. ???. The large item to the left is a frame in which the current conducting wire (1) is suspended. The current from the power supply (2) flows through the centre conductor from the right hand side plate to the left hand side plate, and is returned through the 4 brass bars holding the frame together. This way the contribution of the return current on the measurement is minimised.


Figure 1: *Overall view of the experiment: 1) Current carrying central conductor ("Plasma current"), 2) Power Supply, 3) Single probe, 4) amplifier and integrator, 5) Oscilloscope, 6) Set of 8 pick-up coils, 7) Amplifiers for pick-up coils, 8) DAQ-system ("TU/e nano-Giants"), 9) PC running LabView©.*

A single pick-up coil (3) is shown and a set of 8 pick-up coils mounted onto the central plate. The single probe is connected to an amplifier and analog integrator in a single case (4) and displayed on the scope. The probe can be used to assess the field in the vicinity of the experiment and performance of the analog integrator may be assessed. The set of 8 coils are conditioned by amplifiers (7) and fed into a DAQ system (8) which is controlled by LabView©.


## **2 Detailed description of the experiment**

The aim of this document is to enable other members of the FuseNet Association to duplicate the experiment. The focus of the description is therefore on those items that require own manufacture or modification. The description of the experiment is divided into the following sections:

1. Frame
2. Pick-up coils
3. Current source and connections
4. Amplifiers and integrators


## 2.1 Frame

An overall assembly of the frame is shown in Fig. ??.

Assembly\_Magnetic\_Diagnostic

Figure 2: *Overall assembly*


ALL DIMENSIONS IN MM  
3D CAD MASTER PART NAME:


		Designer: Jeroen BaIm	
		TITLE Assembly: Magnetic Diagnostic Frame	
FIRST ISSUED	DRAWN BY	CHECKED BY	APPROVED BY
SIZE DWG NO. A4		assembly_m_d_drafting	
SCALE 1:1		SHEET 1 OF 1	

1 2 3 4 5 6


ALL DIMENSIONS IN MM  
3D CAD MASTER PART NAME:


Designer : Jeroen Balm

FIRST ISSUED	TITLE	plastic 1 part
DRAWN BY	(File > Properties > Attributes)	
CHECKED BY	SIZE DWG NO.	A4 plasticplate_1
APPROVED BY	SCALE 1:1	SHEET 1 OF 1
	REV.	<input checked="" type="checkbox"/>

ALL DIMENSIONS IN MM  
3D CAD MASTER PART NAME:

 		Designer: Jeroen Balm	
FIRST ISSUED		TITLE	4 parts copper
DRAWN BY			copper
CHECKED BY			(File > Properties > Attributes)
APPROVED BY		SIZE/DWG NO.	A4 copper_bar
		SCALE 1:5	REV. <input checked="" type="checkbox"/>
			SHEET 1 OF 1

10,0

$\phi 10,0$

$\phi 10,0$

276,4

500,0

400,0

400,0

500,0

slot depth 10 mm

Designer : Jeroen Balm

FIRST ISSUED	TITLE
DRAWN BY	Plastic part
CHECKED BY	(File > Properties > Attributes)
APPROVED BY	SIZE DWG NO.
	A4 plasticplate_2
	SCALE 1:1
	SHEET 1 OF 1

REV.

ALL DIMENSIONS IN MM  
3D CAD MASTER PART NAME:

6

5

4

3

2

1

D

C

B

A

D

C

B

A

6


5

4

3



2

1

ALL DIMENSIONS IN MM  
3D CAD MASTER PART NAME:


 		Designer : Jeroen Balm	
FIRST ISSUED		TITLE	
DRAWN BY			copper
CHECKED BY			1 part
APPROVED BY			(File > Properties > Attributes)
		SIZE DWG NO.	copperplate
		A4	
		SCALE 1:1	SHEET 1 OF 1

THIS DRAWING HAS BEEN PRODUCED USING AN EXAMPLE TEMPLATE PROVIDED BY UGS.

FIRST ISSUED	TITLE		
DRAWN BY	Missing <input checked="" type="checkbox"/>		
CHECKED BY	4 parts (File > Properties > Attributes)		
APPROVED BY	SIZE/DWG NO.	thread_m10_parts	REV. <input checked="" type="checkbox"/>
	A4		
	SCALE 1:1		SHEET 1 OF 1

ALL DIMENSIONS IN MM  
3D CAD MASTER PART NAME:

ALL DIMENSIONS IN MM  
3D CAD MASTER PART NAME:

  The PLM Company		Designer: Jeroen Balm	
FIRST ISSUED		TITLE	
DRAWN BY		copper litz	
CHECKED BY		1 part	
APPROVED BY		SIZE DWG NO.	REV.
		A4	copper_litz <input checked="" type="checkbox"/>
		SCALE 1:1	SHEET 1 OF 1

1 2 3 4 5 6

## 2.2 Pick-up coils

One separate pick-up coil is available to freely assess the fields in the vicinity of the experiment while a set of 8 pick-up coils is mounted on the central plate in the frame. Fig. ?? shows the single probe, but mounted on a sliding table, and Fig. ?? shows the set of 8 coils. The engineering drawings are shown on the next pages, excluding the sliding table as this is an undocumented accessory.


Figure 3: *Separate pick-up coil mounted on sliding table.*




Figure 4: *The set of 8 pick-up coils mounted on the central plate in the frame. In the centre the copper litze cable, purposely slightly offset to demonstrate displacement of the current in the plain of the coils.*


Tie rap

 <small>3M Company</small>		THIS DRAWING HAS BEEN PRODUCED USING AN EXAMPLE TEMPLATE PROVIDED BY UGS.	
		TITLE	
FIRST ISSUED	DRAWN BY <i>coilholder</i> <i>assy_coil_coilholder</i>		
CHECKED BY	(File > Properties > Attributes)		
APPROVED BY	SIZE	DWG NO.	REV.
	A4	assy_coil_coilholder	<input checked="" type="checkbox"/>
SCALE 1:1		SHEET 2 OF 2	



ALL DIMENSIONS IN MM  
3D CAD MASTER PART NAME:

SHEET 2 OF 2



Coil:  
 Windings= 1000  
 Wire diameter= 0.15 mm

Designer: Jeroen Balm

FIRST ISSUED	TITLE		
DRAWN BY	Material: PMMA		
CHECKED BY	8 parts		
APPROVED BY	SIZE DWG NO.	REV.	
	A4	spoel_8_parts	<input checked="" type="checkbox"/>
	SCALE 1:1	SHEET 1 OF 1	

ALL DIMENSIONS IN MM  
 3D CAD MASTER PART NAME:


Designer: Jeroen Balm

FIRST ISSUED		TITLE	Material: MS	REV.	<input checked="" type="checkbox"/>
DRAWN BY		8 parts			
CHECKED BY		SIZE DWG NO.	A4	spoolhouder_8_parts	
APPROVED BY		SCALE 1:1		SHEET 1 OF 1	


ALL DIMENSIONS IN MM  
3D CAD MASTER PART NAME:

1 2 3 4 5 6

ALL DIMENSIONS IN MM  
3D CAD MASTER PART NAME:

  The PLM Company		THIS DRAWING HAS BEEN PRODUCED USING AN EXAMPLE TEMPLATE PROVIDED BY UGS.	
FIRST ISSUED	TITLE		
DRAWN BY	Hel Hermann  tyton.com		
CHECKED BY	(File > Properties > Attributes)		
APPROVED BY	SIZE/DWG NO.	REV.	<input checked="" type="checkbox"/>
	A4	grommet	
	SCALE 5:1	SHEET 1 OF 1	

ALL DIMENSIONS IN MM  
3D CAD MASTER PART NAME:

  The PLM Company		THIS DRAWING HAS BEEN PRODUCED USING AN EXAMPLE TEMPLATE PROVIDED BY UGS.	
FIRST ISSUED		TITLE	
DRAWN BY		<input checked="" type="checkbox"/> HellermannTyton.com <small>(File &gt; Properties &gt; Attributes)</small>	
CHECKED BY		SIZE	DWG NO.
APPROVED BY		A4	kabel sleeve
		SCALE 2:1	REV. <input checked="" type="checkbox"/>
			SHEET 1 OF 1

1 2 3 4 5 6      A B C D

### **2.3 Current source and connections**

The power supply has been commercially bought, EA-PSI 9080-100. The unit may be programmed to provide e.g. the shape of a typical plasma pulse. Details of the power supply are attached. The unit was bought by Farnell ([www.farnell.com](http://www.farnell.com)). The connections between the power supply and the frame have been made with flexible copper cable 25 mm<sup>2</sup>. See datasheet attached.

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
- ▶▶ **Output Power:** 1,5kW, 3kW, 6kW and 9kW
- ▶▶ **PFC power factor >0,99**
- ▶▶ **Output Voltages:** 0...80V; 0...300V; 0...600V; 0...750V
- ▶▶ **Remote Sense**
- ▶▶ **Overvoltage Protection OVP**
- ▶▶ **Supervision of  $U_{>(t)}$ ,  $U_{<(t)}$ ,  $I_{>(t)}$  and  $I_{<(t)}$**
- ▶▶ **Alarm Management**
- ▶▶ **Integrated Funktion Generator**
- ▶▶ **User Profile w. Password protection**
- ▶▶ **Block function of the control panel**

### Programmable Power Supplies

The Lab-power series EA-PSI 9000 is based on the series EA-PS 9000. The new conception of the power stage provides to the user either a high output voltage and a high output current. The output values are limited by the output power.

These power supplies are designed for engineering laboratory and complex industrial applications.

The user-friendly and interactive menu guidance simplifies the configuration of user profiles and function sequences. Complex applications can be put into practice very fast and simple.

Function sequences or user profiles can be filed and documented for a better reproducibility of tests.

The integrated monitoring functions with adjustable delays simplify the test assembling because no external monitoring circuits are required.

### Options

- ▶▶ Two Quadrant Module
- ▶▶ Internal resistance regulation
- ▶▶ Communication interface
- ▶▶ CAN BUS interface with gateway-function and adjustable address range
- ▶▶ USB
- ▶▶ RS232
- ▶▶ GPIB (on request)
- ▶▶ Analogue interface, galvanic isolation, with programmable in- and outputs
- ▶▶ Galvanic isolation of the interface inputs and the output for 2kV
- ▶▶ Support of the communication interfaces through LABview VI's
- ▶▶ **EAsyPower**, with the system-, control- and monitoring software **EAsyPower** up to 30 PSI-power supplies can be controlled and monitored at the same time.

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### Display and Control Panel EA-PSI 9000

The display of the power supply shows the:

- Set value of V, A and P
- Actual values of V, A and P
- Operating mode CV, CC or CP
- Set value of internal resistance (optional)
- Output status **ON OFF**
- Access to the power supply
  - local = Front panel operation
  - external = Analogue interface
  - remote = Access via BUS system
- Current Function of the push-buttons
- Alerts, warnings and simple alerts

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
The functions indicated on the display can be activated by actuation of the related push button below the display.

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### User Profiles

Profiles designed to user specifications can be used to set and store user specified parameters. Up to four different profiles can be deposited in the units. Using the system, control and monitoring software **EAsyPower**, the user profiles can be readout, processed and overwritten.

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### Alarm Management

Up to three alarms, warnings or simple alerts can be indicated on the display. Alarms and Warnings can be audibly signaled. Alarms will shut down the output instantly. Warnings must be acknowledged but will not shut down the output. Simple alerts are displayed only.

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All alerts can read out via the optional galvanically isolated analogue interface.

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
### Function Sequence

The function sequence can be controlled via the control panel as well as via the communication interfaces.

The actual values can be read out during a function sequence. They can be stored and later analysed with the system control and monitoring software EAsyPower.

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
### Function characteristic

With the integrated function generator user specified characteristics can be defined according to the function  $f(V, A, \Delta t)$ .

Such function sequence can be repeated as often as wanted.

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A function sequence can be combined of max. five sequences.

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

A sequence is defined via 10 sequence points.

A sequence point is combined of the targeted final value of  $V$  and  $A$  after the time  $\Delta t$ .

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Function characteristics can be read out and stored with the System, control and monitoring software EAsyPower.

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Highly flexible stranded wire with reinforced insulation.

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**Typical Application**

Current feeds and earth/ground wires in machine, plant and accumulator construction. Safety test leads carrying high currents.

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1) UL recognized: File E120880, AWM 3525  
Use: Test Probe Lead up to +60°C


## 2.4 Amplifiers and Integrators

The next pages show in the following order:

- Amplifier and integrator for separate pickup coil
- Amplifiers for set of 8 pickup coils (feeding into the DAQ system)

The signals from the amplifiers of the set of 8 coils are input to a DAQ system such that the signals can be manipulated digitally. We used a locally developed DAQ system for Experiment Automation ("nano-Giant": <http://www.tuedacs.nl/>) but other systems may be used as well.

11 current pulse, XCP1 conversion factor is set to 1000. The 10<sup>3</sup> turns magnetic pick-up coil and the high current wire are represented by a transformer with a poor coupling factor and a 1 to 10<sup>6</sup> inductance ratio. The high current wire has an estimated inductance of 1 nH so the pick-up coil has an inductance of 1 mH. the coupling factor can be adjusted to mimic pick-up coil distance to the main conductor.

The passive low pass input filter is designed for 3.4 kHz

Best common mode suppression by U1 occurs when the ratio  $(R1 + R2)/R5$  match  $(R3 + R4)/R6$ .

A good match can be achieved by manually selecting from batch of 100 pcs of 1% tolerance:

$R1/R3 - 1 < +/- 0.1\%$ ,  $R2/R4 - 1 < +/- 0.1\%$  and  $R5/R6 - 1 < +/- 0.1\%$

The amplification of this stage:  $A(U1) = -R5/(R1 + R2) = -1000$

Integrator stage with U2

$A(U2) = -1/(R7 * C2) = -10000$


A small offset has to be nulled with Potentiometer R9 connected to VDD

VDD -8V or VSS and added to U2+IN through R10

R9: For this simulation the value of R10 is exactly determined (1.763014 M) so that with 0.1 A current in I1 and R9 is at 50% the dnf of the simulation is as low as possible. For the real integrator the value of R10 = 2.2 MOhm

# MagneticInductionIntegrator



Printing Time: April 21, 2013 08:48:35


Traces: Channel A, Channel B

# MagneticInductionIntegrator AC Analysis

Printing Time: April 23, 2013 08:50:17


For input side see MagnIndIntSch text.  
 PreAmp(U1) = -10X low pass 72 Hz  
 PreAmp(U2) = -100X low pass 72 Hz  
 Amp(U3) = 2X low pass 400 Hz  
 The real circuit has +/- 8V stabilisers (78L08 and 79L08). on board, getting power from stabilised +/-1.2V benchtop PS. The 72 Hz low pass filters should be adjusted to slightly under half the sample frequency or the sample frequency should now be set to at least 200 Hz to avoid aliasing.


The setting of potentiometer R10 has no effect in the simulation. The simulation does have an offset from U1 an U2 that I compensated by adding and trimming R12, but R12 is not part of the real circuit.

The choice of OPAGO4AP for U3 was because it is a Hi-Fi Headphone opamp that can drive up to 600 Ohm Load, unlike the LTC1050, allowing short-circuiting the output indefinitely. It also has low distortion and low drift for a non chopper type of opamp.

# MagnIndAmpSch

Printing Time: April 23, 2025 05:00



Traces: Channel A, Channel B

# MagnIndAmpSch AC Analysis

Printing Time: April 23, 2013 15:07:12

### **3 Acknowledgementst**

In our group the experiment has been built by our technical department, Jeroen Balm, Herman de Jong, Gerard Harkema, Jolanda van de Ven, Frans van Setten, and Hans Oosterbeek. We wish to thank the mechanical workshop (GTD), and in particular Mr. Edwin van den Eijnden, for constructive discussions, drawings, and manufacture of parts. The experiment has been built with financial support from the FuseNet Association.