


PRISMAPRO

Communication Possibilities

Agenda

- 1. [Introduction to PrismaPro](#)
- 2. [Functionality](#)
- 3. [Communication possibilities](#)
- 4. [API-Communication](#)
- 5. [EULA](#)
- 6. [PV-Cloud](#)
- 7. [Overview Application](#)
- 8. [Application \(MPP-S7, LabView, ...\)](#)
- 9. [Summary](#)


1. Introduction to PrismaPro



Standard: 1x AI, 2x DI, 1 Relay

Extended: 8x AI, 5x DI, 10x DO, 8x AO
1x total Pressure-Transmitter

2. Functionality



The PrismaPro is a quadrupole mass spectrometer.

The neutral gas particles are ionized in an ion source and the ions are separated from each other in a high-frequency electric quadrupole field according to their mass and charge conditions.

The detection of the filtered ions is carried out using Farady Cup or secondary electron multiplier.


The detected ion flow is proportional to the partial pressure of the respective gas components

We are measuring Ion-Current for adjusted amu


<https://www.youtube.com/watch?v=TrLMNRyQiaA>

3. Communication possibilities

- WebUI:
<http://192.168.1.100>



- PVMassSpec
[Start exe](#)
[YouTube-Video:](#)



3. Communication possibilities

- **MMSP-API:** (Modular Mass Spec - Web-API, http port: 80)
Is a web server, whose API is accessed over HTTP.

API: Application Programming Interface

Is a set of subroutine definitions, protocols and tools for building application software. In general terms, it is a set of clearly defined methods of communication between various software components (*wikipedia*)



[Documentation](#) for the API usually is provided to facilitate usage and Implementation.

https://en.wikipedia.org/wiki/Application_programming_interface

4. API-Communication

- Web-Browser-Input:

<http://192.168.1.100/mmsp/Communication/ipAddress/get>



- cURL (.bat-File)

cURL is a command line tool for getting or sending files using URL syntax (*wikipedia*)

[API-Get IP](#)

[API-Set Sweep .txt](#)

- Access/Excel mit VBA

[API-Access](#)

[API-Excel](#)

4. API-Communication

- JSON (**J**ava**S**cript **O**bject **N**otation)

www.json.org



is an open-standard file format that uses human-readable text to transmit data objects consisting of attribute–value pairs and array data types (*wikipedia*)

It is a syntax-analysis necessary (Parsing) to resolve the dataset into his related elements.

Parse JSON:

Raw-Data: `{"data":"192.168.1.100","name":"got","origin":"/mmsp/communication/ipAddress"}`



data	192.168.1.100
name	got
origin	/mmsp/communication/ipAddress


Example 1: MS-Excel:



PrismaPro_API_Service.xlsm

API-Reference "PrismaPro":		IP-Address PrismaPro: <u>http://192.168.1.100</u>		
<i>The usage of the API is at your own risk !!</i>				API-Reference_EULA.pdf
Nbr	Target	index	set Value:	Title
1	<u>mmssp</u>	get -	-	Modular Mass Spec
2	<u>electronicsInfo</u>	get -		Electronics Information
3	serialNumber	get -	STMS16064	Serial Number
4	monitorFWVersion	get -	-	Monitor Firmware Version
5	controlFWVersion	get -	-	Control Firmware Version
6	genus	get -	4	Genus
7	massRange	get -	-	Mass Range
8	<u>communication</u>	get -	-	Communications Port Configuration
9	ipAddress	get -	192.168.1.100	IP Address
10	macAddress	get -	:A0:41:02:10:	MMSP MAC Address
11	gateway	get -	0.0.0.0	MMSP Gateway
12	mask	get -	255.255.255.0	MMSP Mask

Example 2: MS-Access+VBA



```
Set httpObject = CreateObject("MSXML2.XMLHTTP")
```

```
httpObject.Open "GET", "" & "http://" & IPADDRESS & "/mmsp/communication/ipAddress/get", False
httpObject.Send
```

```
result = httpObject.ResponseText
```

```
MsgBox ("Result of API-Get: " & result)
```



```
Dim JSON As Object
```

```
Set JSON = JsonConverter.ParseJson(result)
```

```
Text3.Value = JSON("data")
```

```
Text6.Value = JSON("name")
```



```
Text8.Value = JSON("origin")
```



Example 3: MS-Excel+VBA:

(Microsoft Excel sheet with Macros)

PrismaPro_API_Service.xlsm



Example 3:

Toggle Emission

```

cmdEmission
Private Sub cmdEmission_Click()
Dim httpObject As Object
Dim i As Integer
Dim IPADDRESS As String
Dim result As String
Dim EMISSION As String

On Error GoTo fault
Set httpObject = CreateObject("MSXML2.XMLHTTP")
Dim JSON As Object

IPADDRESS = Worksheets("API-Ref").Cells(2, 12)

' Emission-----
httpObject.Open "GET", "" & IPADDRESS & "/mmsp/generalControl/setEmission/get", False
httpObject.Send
result = httpObject.ResponseText
Set JSON = JsonConverter.ParseJson(result)
EMISSION = JSON("data")
If (EMISSION = "On") Then
    httpObject.Open "GET", "" & IPADDRESS & "/mmsp/generalControl/setEmission/set?Off", False
    httpObject.Send
Else
    httpObject.Open "GET", "" & IPADDRESS & "/mmsp/generalControl/setEmission/set?On", False
    httpObject.Send
End If

Exit Sub
'Fault????
fault:
    MsgBox ("No/Wrong Answer from PrismaPro! =>Error: " & ERROR)

End Sub
    
```

5. EULA (End User License Agreement)


- Software: API (*limited to ~350*), shall be treated as confidential, ...
- Hardware: PrismaPro
- Advantage for EULA-Users:
API-Reference-List, Login name/password
Examples, Support, free of charge
Free access to PV-Cloud-Folder „EULA-Customer“
- Advantage for Pfeiffer Vacuum :
Better control who uses the API + PrismaPro
Legal certainty, ...

6. PV-Cloud

- Folder: „EULA-Customer“
EULA-Users are invited (by E-Mail)
Password necessary
- Content :
API Reference List
MS-Excel-Sheets with VBA
cURL
API-Batch
LabView-Driver
Flowcharts, ...

Folder will be updated from time to time

6. PV-Cloud




Herunterladen 3 MB

Name	Icon	Größe	Datum
1 Datei	📁	3 MB	vor 2 Monaten
API Reference_EULA.pdf	📄	3 MB	vor 2 Monaten
API-Batch.jpg	🖼️	265 KB	vor 5 Monaten
API-Get_IPAdr.bat	📄	< 1 KB	vor 5 Monaten
Application_API_Description.pdf	📄	710 KB	vor 2 Monaten
Control_PrismaPro.pdf	📄	585 KB	vor 2 Monaten
curl-7.54.1-win64-mingw.zip	📁	2.7 MB	vor 5 Monaten
Firefox_Default_Viewer.png	🖼️	123 KB	vor 2 Monaten
Flowchart_PLC_Communication_ENG.pdf	📄	393 KB	vor 2 Monaten
General_LogIn.png	🖼️	3 KB	vor 5 Monaten


6. PV-Cloud

- Folder: „PrismaPro“
PrismaPro-User are invited
Password necessary
- Content: *Current Software PVMassSpec*
Current Firmware for PrismaPro




Folder will be updated from time to time

7. Overview Applications

- 
■
MPP-S7
PLC: SIEMENS S7-300/400/1500
PC: SIMATIC IPC227E Nanobox PC
Windows 7 / Windows 10

SIEMENS

- 
■
MPP-AB
MPP-IQ-R
MPP-OPC-UA
Rockwell,
Mitsubishi
Beckhoff, Mitsubishi, Rockwell,
(any PLC/PC who can use „OPC-UA-Server“)

**Rockwell
Automation**

 *Allen-Bradley*

 **MITSUBISHI
ELECTRIC**
Changes for the Better

- 
■
LabView-Driver
Free of charge for EULA-Users



8. Application: MPP-S7

Network example 1:

Provide by Pfeiffer Vacuum



192.168.1.100
(Standard-IP-Addr.)


Siemens Box-PC
+ MPP-S7-Application



X1: 192.168.1.3
LAN-2

Customer part

Siemens PLC
S7-300/400/1500



X2: 192.168.1.1




192.168.1.101

(MPP-S7 : Application for Communication Mass-Spectrometer PrismaPro <=> Siemens SIMATIC-S7 PLC)

8. Application: MPP-S7

Network example 2:

SIEMENS



8. Application: MPP-S7

[MPP-S7 Description](#)
[Interface Description](#)

SIEMENS

Main topics:

- Establish a connection to the PrismaPro and to the PLC
- Only work together with PrismaPro! and Siemens-PLC
- Switch On/Off Emission + EM
- Read/Write a recipe from/to the PrismaPro (currently 20 Channels/Recipe)
- Start/Stop a RUN with „Lock“ and „Release“ the PrismaPro for other Applications
- Read the measured values after each RUN (max. 3001 Datapoints/RUN)
- Cyclic reading: Internal + external Pressure (mBar & Pas & Torr)
Faults and Warnings
State of Emission, EM, Scan,
Digital/Analog Inputs + Outputs (Option)

Feature: Adjust a recipe in PVMassSpec and use it in the Application

8. Application: MPP-S7

PLC 1516-3 PN/DP PLC Online (Watch-List):



Name	Adresse	Anzeigeformat	Beobachtungswert	Steu...	...	Variablen-Kommentar
"PrismaPro".Common.Lifebit	%DB100.DBX0.0	BOOL	<input checked="" type="checkbox"/> TRUE	TRUE	<input type="checkbox"/>	Lifebit, is set all 1s by PLC, Reset by PrismaPro
"PrismaPro".Common.PrismaProErrorCode	%DB100.DBD2	DEZ	0	0	<input type="checkbox"/>	Errorcode from PrismaPro Device
"PrismaPro".Common.DriverErrorCode	%DB100.DBD6	DEZ	0	0	<input type="checkbox"/>	Errorcode from PrismaPro Device Driver
"PrismaPro".DeviceInfo.Genus	%DB100.DBB10	DEZ	4		<input type="checkbox"/>	PrismaPro: Communication Device Type: 4=PrismaPro
"PrismaPro".DeviceInfo.massRange	%DB100.DBW12	DEZ	100		<input type="checkbox"/>	PrismaPro: Possible mass Range
"PrismaPro".DeviceInfo.SensorName	P#DB100.DBX1...	String	'QMSPrismaPro'		<input type="checkbox"/>	PrismaPro: Sensor Name
"PrismaPro".Communication.IPAdress	P#DB100.DBX3...	String	'172.16.14.42'	''	<input type="checkbox"/>	PrismaPro: PrismaPro Device IP Adress
"PrismaPro".Communication.Control	%DB100.DBW54	DEZ	256		<input type="checkbox"/>	PrismaPro: controlled by PrismaPro Driver; 0=none,1=self,2=other
"PrismaPro".Control.SetEmissionOn	%DB100.DBX58.1	BOOL	<input type="checkbox"/> FALSE	TRUE	<input type="checkbox"/>	Set emission State On (PLC 0->1, 1->0 by Driver)
"PrismaPro".Control.SetEmissionOff	%DB100.DBX58.2	BOOL	<input type="checkbox"/> FALSE	TRUE	<input type="checkbox"/>	Set emission State Off (PLC 0->1, 1->0 by Driver)
"PrismaPro".Control.SetEMOn	%DB100.DBX58.3	BOOL	<input type="checkbox"/> FALSE	TRUE	<input type="checkbox"/>	Set electron Multiplier State On (PLC 0->1, 1->0 by Driver)
"PrismaPro".Control.SetEMOff	%DB100.DBX58.4	BOOL	<input type="checkbox"/> FALSE	TRUE	<input type="checkbox"/>	Set electron Multiplier State Off (PLC 0->1, 1->0 by Driver)
"PrismaPro".Control.StartScan	%DB100.DBX58.5	BOOL	<input type="checkbox"/> FALSE	TRUE	<input checked="" type="checkbox"/>	Start Measurement (PLC 0->1, 1->0 by Driver)
"PrismaPro".Control.StopScan	%DB100.DBX58.6	BOOL	<input type="checkbox"/> FALSE	TRUE	<input type="checkbox"/>	Stop Measurement (PLC 0->1, 1->0 by Driver)
"PrismaPro".Control.SaveSetup	%DB100.DBX58.7	BOOL	<input type="checkbox"/> FALSE	TRUE	<input type="checkbox"/>	Save Scan Setup (PLC 0->1, 1->0 by Driver) => PLC to PrismaPro
"PrismaPro".Control.ReadHardware	%DB100.DBX58.8	BOOL	<input type="checkbox"/> FALSE	FALSE	<input type="checkbox"/>	Switch for Hardware IO reading (0 = off, 1 = on)

8. Application: MPP-S7

PLC 1516-3 PN/DP PLC Online (DB100 Part 1)



PrismaPro					
	Name	Datentyp	Offset	Beobachtungswert	Kommentar
1	Static				
2	Common	Struct	0.0		Common Control and State
3	Lifebit	Bool	0.0	TRUE	Lifebit, is set all 1s by PLC, Reset by PrismaPro
4	PrismaProErrorCode	DInt	2.0	0	Errorcode from PrismaPro Device
5	DriverErrorCode	DInt	6.0	0	Errorcode from PrismaPro Device Driver
6	DeviceInfo	Struct	10.0		Electronic Info from Prisma Device
7	Genus	Byte	10.0	16#04	PrismaPro: Communication Device Type: 4=PrismaPro
8	massRange	Word	12.0	16#0064	PrismaPro: Possible mass Range
9	SensorName	String[20]	14.0	'QMSPrismaPro'	PrismaPro: Sensor Name
10	Communication	Struct	36.0		Communication Info from Prisma Device
11	IPAdress	String[15]	36.0	'172.16.14.42'	PrismaPro: PrismaPro Device IP Adress
12	Control	Int	54.0	256	PrismaPro: controlled by PrismaPro Driver; 0=none,1=self,2=other
13	ControllInfo	Int	56.0	0	PrismaPro: controlled by external device; not used
14	Control	Struct	58.0		=== PrismaPro Control ===
15	SetEmissionOn	Bool	58.0	FALSE	Set emission State On (PLC 0->1, 1->0 by Driver)
16	SetEmissionOff	Bool	58.1	FALSE	Set emission State Off (PLC 0->1, 1->0 by Driver)
17	SetEMOn	Bool	58.2	FALSE	Set electron Multiplier State On (PLC 0->1, 1->0 by Driver)
18	SetEMOff	Bool	58.3	FALSE	Set electron Multiplier State Off (PLC 0->1, 1->0 by Driver)
19	StartScan	Bool	58.4	FALSE	Start Measurement (PLC 0->1, 1->0 by Driver)
20	StopScan	Bool	58.5	FALSE	Stop Measurement (PLC 0->1, 1->0 by Driver)
21	SaveSetup	Bool	58.6	FALSE	Save Scan Setup (PLC 0->1, 1->0 by Driver) => PLC to PrismaPro
22	ReadHardware	Bool	58.7	FALSE	Switch for Hardware IO reading (0 = off, 1 = on)
23	UserTimer1	Time	60.0	T#0MS	Set User Timer 1 (Set value by PLC)
24	UserTimer2	Time	64.0	T#0MS	Set User Timer 2 (Set value by PLC)

8. Application: MPP-S7

PLC 1516-3 PN/DP PLC Online (DB100 Part 2)



[-] ▶ In_ScanSetup	Struct	98.0		=== Scan Setup === PLC to PrismaPro
[-] ▶ Out_ScanSetup	Struct	6920.0		=== Feedback of Scan Setup === PrismaPro to PLC
[-] ▶ IOState	Struct	13742.0		PrismaPro: State of the iterated IO signals
[-] ▾ Scans	Struct	13838.0		=== Single Channel Scan results ===
[-] ▫ DateTime	Date_And_Time	13838.0	DT#2019-03-11-09:50:06.10	PrismaPro: Date and Time of last Scan - Current DateTime of PC (UTC)
[-] ▫ ▾ Result	Array[0..3000] of Struct	13846.0		PrismaPro: Measured Channel Value
[-] ▫ ▫ ▾ Result[0]	Struct	13846.0		PrismaPro: Measured Channel Value
[-] ▫ ▫ ▫ Mass	Real	13846.0	0.0	PrismaPro: Measured Mass
[-] ▫ ▫ ▫ Value	Real	13850.0	4.760849E-12	PrismaPro: Measured Mass Value
[-] ▫ ▫ ▾ Result[1]	Struct	13854.0		PrismaPro: Measured Channel Value
[-] ▫ ▫ ▫ Mass	Real	13854.0	0.2	PrismaPro: Measured Mass
[-] ▫ ▫ ▫ Value	Real	13858.0	2.154627E-12	PrismaPro: Measured Mass Value
[-] ▫ ▫ ▾ Result[2]	Struct	13862.0		PrismaPro: Measured Channel Value
[-] ▫ ▫ ▫ Mass	Real	13862.0	0.4	PrismaPro: Measured Mass
[-] ▫ ▫ ▫ Value	Real	13866.0	1.589171E-12	PrismaPro: Measured Mass Value
[-] ▫ ▫ ▾ Result[3]	Struct	13870.0		PrismaPro: Measured Channel Value
[-] ▫ ▫ ▫ Mass	Real	13870.0	0.6	PrismaPro: Measured Mass
[-] ▫ ▫ ▫ Value	Real	13874.0	1.703431E-12	PrismaPro: Measured Mass Value
[-] ▫ ▫ ▾ Result[4]	Struct	13878.0		PrismaPro: Measured Channel Value
[-] ▫ ▫ ▫ Mass	Real	13878.0	0.8	PrismaPro: Measured Mass
[-] ▫ ▫ ▫ Value	Real	13882.0	2.200243E-12	PrismaPro: Measured Mass Value
[-] ▫ ▫ ▶ Result[5]	Struct	13886.0		PrismaPro: Measured Channel Value
[-] ▫ ▫ ▶ Result[6]	Struct	13894.0		PrismaPro: Measured Channel Value




8. Application: MPP-OPC-UA




Network concept:

(wikipedia: [OpenPlatformCommunications-UnifiedArchitecture](https://en.wikipedia.org/wiki/OpenPlatformCommunications-UnifiedArchitecture))



8. Application: MPP-IQ-R

Network concept study:





8. Application: LabView-Driver+Example

Frontpanel von PV PrismaPro Single Measurement.vi auf Kalibrierpumpstand.lvproj/Mein Computer

Device Configuration | Scan Setup | Scans

IP address
192.168.1.100
username
OEM-33-API
password
8c38bf2d-f6db2426-ca6aa865-bb0c8424

connected

Connect

Emission

Electron Multiplier

Stop Scan

Linear Scale

len Current [A]

Sample

1E-7
1E-6
1E-5
1E-4
1E-3
1E-2
1E-1
1E-0
1E-1
1E-2
1E-3
1E-4
1E-5
1E-6
1E-7
1E-8
1E-9
1E-10
1E-11
1E-12
1E-13
1E-14
1E-15

0 5 10 15 20 25 30 35 40 45 50 55 60 65 70 75 80 85 90 95 100

1. Tab Device Configuration
Set the proper IP address of the connected PrismaPro.
If necessary modify login-information (username/password).
Press the "Connect" button.

2. Tab Scan Setup
Modify to the desired scan setup and send the final values to the PrismaPro by pressing the "Write ScanSetup" button.
Note that you have to modify the stopChannel if you configure additional channel.
Setting the scanCount to -1 tells the instrument to continue scanning indefinitely.

3. Tab Scans
Turn on the emission.
Note that setting the emission on can take several seconds.
WARNING:
The sensor **MUST** be under vacuum when emission is turned on, or the filament will be destroyed.
Turn on the electron multiplier.
Start the scan-process by pressing the "Start Scan" button.

9. Summary

- For this Communication, a **PrismaPro** is required.
- For general Application and Calibration use **PVMassSpec**
- For a higher-level Communication (e.g. LabView, PC, ...) use additional the **MMSP-API**
- For the Communication with the API, the signed **EULA** is necessary but not for the MPP-S7-Application.



Thank you very much for your attention!

