TOKAMAK GOLEM FOR FUSION EDUCATION - CHAPTER 11

P. Macha^{1,2}, K. Hromasova^{1,2}, D. Kropackova¹, M. Lauerova³, A. Socha⁴, J. Malinak¹, L. Lobko¹, D. Cipciar⁵, J. Adamek², J. Cecrdle¹, V. Svoboda¹, J. Stockel²

Faculty of Nuclear Sciences and Physical Engineering CTU in Prague, Czech Rep., ² Institute of Plasma Physics AS CR, Prague, Czech Rep., ³ Gymnazium Nad Aleji 1952, Praha, Czech Rep., ⁴ Gymnazium a Strední odborna skola Cihelni 410, Frydek-Mistek, Czech Rep., ⁵ Faculty of Science, Masaryk University, Brno, Czech Rep.


The GOLEM tokamak



- Parameters: $B_{\rm t} < 0.5$ T, $I_{\rm p} < 8$ kA, pulse length < 15 ms.
- An educational device for domestic as well as for foreign students via remote participation/handling [1].
- Students become familiar with probe measurements, data analysis and basic tokamak diagnostics.
- Subject of Bachelor's degree projects and Master's degree theses.

External plasma stabilization


- 2 external windings generating predefined poloidal magnetic field used for plasma control (horizontal, vertical).
- Plasma position is determined by Mirnov coils.
- The influence of the horizontal mg. field on the discharge duration is shown.



Plasma vertical displacement for discharges with stabilization generating horizontal magnetic field.

Long-range correlations

- Long-range correlations studies using 2 toroidaly separated Langmuir probes with focus on GAMs.
- The coherence for GAM confirmation was not reached yet.



Example of one of the long-range coherence in the edge plasma of GOLEM tokamak.

Lithization tests in vacuum tube

- Lithization setup tests performed in a small vacuum tube.
- Several critical problems need to be handled before usage in the tokamak wall - oxidation of the metallic lithium and potential melting of electrodes.
- A spectral line of neutral Li has been observed.




Snapshot of the glow discharge with the apparent red neutral lithium line.

References


- [1] V. Svoboda, et al., Fus. Eng. and Des. 68,
- 1310-1314 (2011) [2] J. Adamek et al 2021 Nucl. Fusion 61 036023 https://iopscience.iop.org/article/10.1088/1741-4326/abd41d
- [3] J. Adamek et al., Czechoslovak Journal of Physics. (2004) 54, 95–99.

Fast ion temperature measurement using swept ball-pen probe


- Ion temperature is measured with 5 μs temporal resolution based on the measurements of the electron branch of a ball-pen probe (BPP) IV characteristics [2].
- The probe collector is biased with a voltage swept between -30V to +130 V at a frequency of 100 kHz.
- The $T_{\rm i}$ is obtained from $I(V) = I_{\rm sat}^+ \cdot \left(\exp(\alpha_{\rm BPP}) \cdot [1 + K \cdot (V \Phi)] \exp\left(\frac{\Phi V}{T_{\rm i}}\right) \right)$, $\alpha = \ln \frac{I_{\rm sat}^-}{I_{\rm i}^+} = 0.25 \pm 0.09 \, (B_{\rm t} > 0.22 \, {\rm T})$.
- Cut-off fitting technique is applied to all the IV characteristics.
- Fluctuations of the ion temperature ranging between 5 eV up to 40 eV reveal the turbulent behavior of the edge plasma.
- NON-Gaussian shaped histograms of T_e and T_e are observed with a peak at low temperature and a tail towards high temperatures.



Cut-of techniqure.



Temporal evolution of T_i and T_e .




Histograms of T_i and T_e .

Electron temperature measurements using rail probe

- The rail probe concept can sustain exceptionally high heat flux and reduce the sheath expansion effect.
- A probe head consists of a rail probe (RP, length = 40 mm, wide = 2 mm), Langmuir probe (LP, length 1.5 mm, diameter $\frac{1}{2}$ 20 -1 mm), and ball-pen probe (BPP) [3] has been designed
- Special manipulator with changable inclination to B_t within $\pm 10^{\circ}$ has been installed.
- Electron temperature is measured using a swept Langmuir and rail probe (f = 5 kHz) and a floating ball-pen probe.
- Capability of RP to reduce the sheath expansion effect was confirmed. Good agreement between LP, RP and BPP electron tempera-

ture measurements for large magnetic field.

Left) Comparison of $T_{\rm e}$ measured by BPP, LP and RP. Right)








Diagram of the combined probe head.

Measurements of HXR radiation

- Scintilation detectors were used for HXR spectromery.
- Electronics of standard photomultiplier tubes of scintillation detectors does not withstand too large HXR fluxes
- Two problems occured:
 - Standard photomultiplier tubes of scintilation detectors can not withstand high HXR (NaI(Tl) detector drops around 8 ms).
 - Pile-up areas and areas of individual peaks.
- Optimal setup must be found by ensuring sufficient lead shielding and the distance from tokamak.




Comparison of HXR signals from 4 different scintillation detectors.



Comparison of piled-up signals and individual peaks.

Turbulent structures


- Exchange turbulence (blob-hole pair generation and propagation) in the plasma edge enhances energy and particle losses.
- Double rake probe (tokamak bottom port) measured I_{sat} at r = 37 - 90 mm (limiter at r = 85 mm).
- \bullet $I_{\rm sat}$ histograms found asymmetric with positive skewness indicates the presence of blobs.
- Skewness seems to decrease to negative values at r = 40 mm, possible location of the blob birth zone.



Radial profile of ion saturated current skewness. Positive values indicate the presence of blobs throughout the investigated region.

Electron temperature measurements

- Swept Langmuir probe verified but complicated and timeconsuming.
- Combined method (floating ball-pen and Langmuir probe) straightforward and high time resolution, but rather new.
- It was verified that both methods give the same results.
- We suggest that the combined method is suitable for measuring the edge plasma $T_{\rm e}$.



Time evolution of $T_{\rm e}$ in two identical GOLEM discharges, showing good correspondence between the two methods.

Acknowledgment

This work was supported by the Grant Agency of the Czech Technical University in Prague, grant No. SGS19/180/OHK4/3T/14.

Contact us

Tokamak GOLEM, Břehová 7, Prague, Czech Republic golem.fjfi.cvut.cz svoboda@fjfi.cvut.cz