

Title

Introduction to the tokamak operation (GOLEM specific) - Level 1

Vojtěch Svoboda
on behalf of the tokamak GOLEM team
for **GOMTRAIC** training session

2019-03-04

Table of Contents

1 Introduction

2 Data handling @ the Tokamak GOLEM

3 Conclusion

4 Appendix

Table of Contents

1 Introduction

- The GOLEM tokamak - introduction
- The (GOLEM) tokamak concept
- The scenario to make the (GOLEM) tokamak discharge
- The scenario to discharge virtually
- The GOLEM tokamak - guide tour
- The GOLEM tokamak - basic diagnostics


2 Data handling @ the Tokamak GOLEM

3 Conclusion

4 Appendix

The GOLEM tokamak basic characteristics

The grandfather of all tokamaks (ITER newsline 06/18)



- Vessel major radius $R_0 = 0.4$ m
- Vessel minor radius $r_0 = 0.1$ m
- Plasma minor radius: $a \approx 0.06$ m
- Toroidal magnetic field: $B_t < 0.5$ T
- Plasma current: $I_p < 8$ kA
- Electron density:
 $n_e \approx 0.2 - 3 \times 10^{19} \text{ m}^{-3}$
- Effective ion charge: $Z_{\text{eff}} \approx 2.5$
- Electron temperature: $T_e < 100$ eV
- Ion temperature: $T_i < 50$ eV
- Discharge duration: $\tau_p < 25$ ms
- (Electron) energy confinement time:
 $\tau_e \approx 50$ us

The GOLEM tokamak for education - historical background

Kurchatov Institute near Moscow,
Soviet Union
1960: **TM1-MH**



1974

Culham Centre for Fusion Energy
Great Britain
1989: **COMPASS-D**



2006

Institute of Plasma Physics
Czech republic
CASTOR **COMPASS**



2008

Czech Technical University Prague
Czech republic
GOLEM



GOLEM

... somewhere, in the ancient cellars of Prague,

there is hidden indeed "infernal" power. Yet it is the very power of celestial stars themselves. Calmly dormant, awaiting mankind to discover the magic key, to use this power for their benefit...



At the end of the 16th century, in the times when the Czech lands were ruled by Emperor Rudolf II, in Prague, there were Rabbi Judah Loew, well known alchemist, thinker, scholar, writer and inventor of the legendary GOLEM - a clay creature inspired with the Universe power that pursued his master's command after being brought to life with a shem, . Golem is not perceived as a symbol of evil, but rather as a symbol of power which might be useful but is very challenging to handle. To learn more of the Golem legend, see e.g. Wikipedia/Golem.

The global schematic overview of the GOLEM experiment




Table of Contents

1 Introduction


- The GOLEM tokamak - introduction
- The (GOLEM) tokamak concept
- The scenario to make the (GOLEM) tokamak discharge
- The scenario to discharge virtually
- The GOLEM tokamak - guide tour
- The GOLEM tokamak - basic diagnostics

2 Data handling @ the Tokamak GOLEM


3 Conclusion

4 Appendix


Our goal: the technology to create a μ Sun on the Earth




Magnetic confinement requires toroidal geometry




A chamber contains the thermonuclear reaction



Toroidal magnetic field coils confine the plasma



A transformer action creates and heats the plasma



The final technology altogether




Table of Contents

1 Introduction

- The GOLEM tokamak - introduction
- The (GOLEM) tokamak concept
- The scenario to make the (GOLEM) tokamak discharge
- The scenario to discharge virtually
- The GOLEM tokamak - guide tour
- The GOLEM tokamak - basic diagnostics


2 Data handling @ the Tokamak GOLEM

3 Conclusion

4 Appendix

Plasma in Tokamak (GOLEM) - the least to do

To do:




- session start phase:
 - Evacuate the chamber
- pre-discharge phase
 - Charge the capacitors
 - Fill in the working gas
 - Preionization
- discharge phase
 - Toroidal magnetic field to confine plasma
 - Toroidal electric field to breakdown neutral gas into plasma
 - Toroidal electric field to heat the plasma
 - Plasma positioning
 - Diagnostics
- post-discharge phase

Plasma in Tokamak (GOLEM) - the least to do

To do:


- session start phase:
 - Evacuate the chamber
- pre-discharge phase
 - Charge the capacitors
 - Fill in the working gas
 - Preionization
- discharge phase
 - Toroidal magnetic field to confine plasma
 - Toroidal electric field to breakdown neutral gas into plasma
 - Toroidal electric field to heat the plasma
 - Plasma positioning
 - Diagnostics
- post-discharge phase



Plasma in Tokamak (GOLEM) - the least to do


To do:

- session start phase:
 - **Evacuate the chamber**
- pre-discharge phase
 - Charge the capacitors
 - Fill in the working gas
 - Preionization
- discharge phase
 - Toroidal magnetic field to confine plasma
 - Toroidal electric field to breakdown neutral gas into plasma
 - Toroidal electric field to heat the plasma
 - Plasma positioning
 - Diagnostics
- post-discharge phase



Plasma in Tokamak (GOLEM) - the least to do


To do:



- session start phase:
 - Evacuate the chamber
- pre-discharge phase
 - **Charge the capacitors**
 - Fill in the working gas
 - Preionization
- discharge phase
 - Toroidal magnetic field to confine plasma
 - Toroidal electric field to breakdown neutral gas into plasma
 - Toroidal electric field to heat the plasma
 - Plasma positioning
 - Diagnostics
- post-discharge phase


Plasma in Tokamak (GOLEM) - the least to do

To do:



- session start phase:
 - Evacuate the chamber
- pre-discharge phase
 - Charge the capacitors
 - **Fill in the working gas**
 - Preionization
- discharge phase
 - Toroidal magnetic field to confine plasma
 - Toroidal electric field to breakdown neutral gas into plasma
 - Toroidal electric field to heat the plasma
 - Plasma positioning
 - Diagnostics
- post-discharge phase

Plasma in Tokamak (GOLEM) - the least to do




To do:

- session start phase:
 - Evacuate the chamber
- pre-discharge phase
 - Charge the capacitors
 - Fill in the working gas
 - **Preionization**
- discharge phase
 - Toroidal magnetic field to confine plasma
 - Toroidal electric field to breakdown neutral gas into plasma
 - Toroidal electric field to heat the plasma
 - Plasma positioning
 - Diagnostics
- post-discharge phase

Plasma in Tokamak (GOLEM) - the least to do


To do:



- session start phase:
 - Evacuate the chamber
- pre-discharge phase
 - Charge the capacitors
 - Fill in the working gas
 - Preionization
- discharge phase
 - **Toroidal magnetic field to confine plasma**
 - Toroidal electric field to breakdown neutral gas into plasma
 - Toroidal electric field to heat the plasma
 - Plasma positioning
 - Diagnostics
- post-discharge phase

Plasma in Tokamak (GOLEM) - the least to do

To do:



- session start phase:
 - Evacuate the chamber
- pre-discharge phase
 - Charge the capacitors
 - Fill in the working gas
 - Preionization
- discharge phase
 - Toroidal magnetic field to confine plasma
 - **Toroidal electric field to breakdown neutral gas into plasma**
 - **Toroidal electric field to heat the plasma**
 - Plasma positioning
 - Diagnostics
- post-discharge phase

Tokamak GOLEM - schematic experimental setup




Table of Contents

1 Introduction


- The GOLEM tokamak - introduction
- The (GOLEM) tokamak concept
- The scenario to make the (GOLEM) tokamak discharge
- The scenario to discharge virtually
- The GOLEM tokamak - guide tour
- The GOLEM tokamak - basic diagnostics

2 Data handling @ the Tokamak GOLEM


3 Conclusion

4 Appendix


Introduce the working gas (Hydrogen x Helium)




Switch on the preionization



Introduce the magnetic field



Introduce the electric field



Plasma ..




Table of Contents

1 Introduction

- The GOLEM tokamak - introduction
- The (GOLEM) tokamak concept
- The scenario to make the (GOLEM) tokamak discharge
- The scenario to discharge virtually
- The GOLEM tokamak - guide tour
- The GOLEM tokamak - basic diagnostics

2 Data handling @ the Tokamak GOLEM


3 Conclusion

4 Appendix


Infrastructure room (below tokamak) 10/16




Infrastructure room (below tokamak) 10/16




Tokamak room (North) 10/16



Tokamak room (North) 10/16



Tokamak room (South) 10/16



Tokamak room (South) 10/16




Table of Contents

1 Introduction


- The GOLEM tokamak - introduction
- The (GOLEM) tokamak concept
- The scenario to make the (GOLEM) tokamak discharge
- The scenario to discharge virtually
- The GOLEM tokamak - guide tour
- The GOLEM tokamak - basic diagnostics

2 Data handling @ the Tokamak GOLEM




3 Conclusion

4 Appendix



The GOLEM tokamak - basic diagnostics






Loop voltage U_l @ the GOLEM tokamak




Toroidal magnetic field B_t @ the tokamak GOLEM



Total current I_{ch+p}



Basic diagnostics traces at the GOLEM tokamak



Remote operation web app - Control room

GOLEM remote Introduction **Control room** Live Results Prague Access: Level 1 Help

Introduction Working gas Preionization Magnetic field Current drive Submit

This web interface will walk you through the process of configuring a discharge in the GOLEM tokamak. All settable values are perfectly safe. Proceed through each step by setting the desired values and then clicking the **Next** button. You can always go to a specific step by clicking its tab.

Preionization (electron gun)

Vacuum stand

Toroidal magnetic field

Current drive

GAS handling H_2/He

2kV

$C_g = 23 \text{ mF}$

$C_0 = 11.3 \text{ mF}$

3D model rendering method: **Static image (fast)** Interactive X3DOM (slower)

Next

Table of Contents

- 1** Introduction
- 2** Data handling @ the Tokamak GOLEM
- 3** Conclusion
- 4** Appendix

GOLEM basic Data Acquisition System (DAS)


- $U_I, U_{B_t}, U_{I_{p+ch}}, I_{rad}$
- $\Delta t = 1\mu s/f = 1MHz$.
- Integration time = 40 ms, thus DAS produces 6 columns x 40000 rows data file.
- Discharge is triggered at 5th milisecond after DAS to have a zero status identification.




Data file example, DAS $\Delta t = 1\mu s/f = 1MHz$ (neutral gas into plasma breakdown focused)

t	$\approx U_I$	$\approx \frac{U_{dB_T}}{dt}$	$\approx \frac{U_d(I_{p+ch})}{dt}$	$\approx I_{rad}$
first	\approx	7405	lines ..	
:	:	:	:	:
0.007383	1.53931	0.390015	0.048828	0.001831
0.007384	1.53686	0.395508	0.067749	0.00061
0.007385	1.54053	0.391235	0.079956	0.00061
0.007386	1.53686	0.38147	0.072632	0
0.007387	1.54297	0.397949	0.059204	0.00061
0.007388	1.54053	0.384521	0.05249	0.00061
0.007389	1.54053	0.39856	0.068359	0.001221
0.00739	1.54053	0.393677	0.082397	0.001221
0.007391	1.53809	0.38208	0.072632	0.001221
0.007392	1.54297	0.400391	0.056763	0.00061
0.007393	1.54419	0.383911	0.053101	0.00061
0.007394	1.53931	0.397339	0.068359	0.001221
0.007395	1.54297	0.391846	0.084229	0.00061
0.007396	1.54541	0.394897	0.074463	0.00061
0.007397	1.54297	0.388184	0.056763	0.001221
0.007398	1.54297	0.391846	0.056763	0.00061
0.007399	1.54297	0.394287	0.06897	0.00061
:	:	:	:	:
next	\approx	32500	lines ..	
:	:	:	:	:
:	:	:	:	:


Plot #4665 U_l graph




python




gnuplot




mathematica



idl



octave



matlab

Data access

All the recorded data and the settings for each discharge (shot) are available at the GOLEM website. The root directory for the files is:

`http://golem.fjfi.cvut.cz/shots/<#ShotNo>/`

The most recent discharge has the web page:

`http://golem.fjfi.cvut.cz/shots/0`

Particular data from DAS or specific diagnostics have the format:

`http://golem.fjfi.cvut.cz/utils/data/<#ShotNo>/<identifier>`

An overview of available data with identifiers, units, description, etc. for each discharge is at

`http://golem.fjfi.cvut.cz/shots/<#ShotNo>/Data.php`

Matlab

```
ShotNo=22471;
baseURL='http://golem.fjfi.cvut.cz/utils/data/';
identifier='loop_voltage';
%Create a path to data
dataURL=strcat(baseURL,int2str(ShotNo), '/', identifier);
% Write data from GOLEM server to a local file
urlwrite(dataURL, identifier);
% Load data
data = load(identifier, '\t');
% Plot and save the graph
plot(data(:,1)*1000, data(:,2), '.');
xlabel('Time [ms]')
ylabel('U_I [V]')
saveas(gcf, 'plot', 'jpeg');
exit;
```

Jupyter (python)

```
import numpy as np
import matplotlib.pyplot as plt

shot_no = 22471
identifier = "loop_voltage"
# create data cache in the 'golem_cache' folder
ds = np.DataSource('golem_cache')
#Create a path to data and download and open the file
base_url = "http://golem.fjfi.cvut.cz/utils/data/"
data_file = ds.open(base_url+str(shot_no)+'/+'+identifier)
#Load data from the file and plot to screen and to disk
data = np.loadtxt(data_file)
plt.plot(data[:,0], data[:,1]) #1. column vs 2. column
plt.savefig('graph.jpg')
plt.show()
```

Gnuplot

```
set macros;
ShotNo = "22471";
baseURL = "http://golem.fjfi.cvut.cz/utils/data/";
identifier = "loop_voltage";
#Create a path to data
DataURL= "@baseURL@ShotNo/@identifier";
#Write data from GOLEM server to a local file
!wget -q @DataURL;
#Plot the graph from a local file
set datafile separator "\t";
plotstyle = "with_lines_linestyle_-1"
plot 'loop_voltage' using 1:2 @plotstyle;
exit;


# command line execution:
# gnuplot Uloop(gp -persist
```

GNU Wget

GNU Wget is a free software package for retrieving files using HTTP, HTTPS and FTP, the most widely-used Internet protocols. It is a non-interactive commandline tool, so it may easily be called from scripts, cron jobs, terminals without X-Windows support, etc.

- Runs on most UNIX-like operating systems as well as Microsoft Windows.
- Homepage: <http://www.gnu.org/software/wget/>
- Basic usage:
 - To get U_i : wget http://golem.fjfi.cvut.cz/utils/data/<\#ShotNo>/loop_voltage
 - To get whole shot: wget -r -nH --cut-dirs=3 --no-parent -l2 -Pshot http://golem.fjfi.cvut.cz/shots/<\#ShotNo>

Excel



File → Open →

<http://golem.fjfi.cvut.cz/utils/data/<#ShotNo>/<identifier>>

Spreadsheets (Excel and others)

are not recommended, only tolerated.

Table of Contents


1 Introduction

2 Data handling @ the Tokamak GOLEM

3 Conclusion

4 Appendix

The global schematic overview of the GOLEM experiment




Production

- Everything via http://golem.fjfi.cvut.cz/GOMTRAIC_
 - This presentation
 - Control rooms
 - Contact: Vojtech Svoboda,
+420 737673903,
svoboda@fjfi.cvut.cz
 - Chat:
tokamak.golem@gmail.com or
skype: tokamak.golem



Recommended values for the GOLEM tokamak operation



- Preionization: Top electron gunn
- Gas: Hydrogen. A Working gas pressure: p_{WG} [mPa] $\in < 0, 40 >$ mPa
- A voltage to charge the Current drive field E_t capacitor: U_{E_t} [V] $\in < 400, 700 >$ V
- A voltage to charge the Toroidal magnetic field B_t capacitor: U_{B_t} [V] $\in < 600, 1200 >$ V
- Time delay of the E_t trigger with respect to the B_t trigger: T_{CD} [μ s] $\in < 0, 10000 >$ μ s

Thank you for your attention

Tokamak TM1

@Kurchatov Institute near Moscow
~1960-1977



SCIENCE

Tokamak CASTOR


@Institute of Plasma Physics, Prague
1977-2007



SCIENCE
& education

Tokamak GOLEM


@Czech Technical University, Prague
2007-



EDUCATION
& science

... with the biggest
control room
in the world ..

Tokamak Golem **REMOTE** for MASTER (Level I)
The smallest & oldest operational tokamak with the biggest control room in the world



Acknowledgement

Financial support highly appreciated:

CTU RVO68407700, SGS 17/138/OHK4/2T/14, GAČR GA18-02482S,
EU funds CZ.02.1.01/0.0/0.0/16_019/0000778 and
CZ.02.2.69/0.0/0.0/16_027/0008465, IAEA F13019, FUSENET and
EUROFUSION.

Students, teachers, technicians (random order):

Vladimír Fuchs, Ondřej Grover, Jindřich Kocman, Tomáš Markovič, Michal Odstrčil, Tomáš Odstrčil, Gergo Pokol, Igor Jex, Gabriel Vondrášek, František Žácek, Lukáš Matěna, Jan Stockel, Jan Mlynář, Jaroslav Krbec, Radan Salomonovič, Vladimír Linhart, Kateřina Jiráková, Ondřej Ficker, Pravesh Dhyani, Juan Ignacio Monge-Colepicolo, Jaroslav Čeřovský, Bořek Leitl, Martin Himmel, Petr Švihra, Petr Mácha, Vojtěch Fišer, Filip Papoušek, Sergei Kulkov, Martin Imríšek.

Table of Contents

1 Introduction

2 Data handling @ the Tokamak GOLEM

3 Conclusion

4 Appendix

References I

-  V. Svoboda, B. Huang, J. Mlynar, G.I. Pokol, J. Stockel, and G Vondrasek. Multi-mode Remote Participation on the GOLEM Tokamak. *Fusion Engineering and Design*, 86(6-8):1310–1314, 2011.
-  Brotankova, J. Study of high temperature plasma in tokamak-like experimental devices. PhD. thesis 2009.
-  Tokamak GOLEM team. Tokamak GOLEM at the Czech Technical University in Prague. <http://golem.fjfi.cvut.cz>, 2007. [Online; accessed March 2, 2019].
-  J. Wesson. *Tokamaks*, volume 118 of *International Series of Monographs on Physics*. Oxford University Press Inc., New York, Third Edition, 2004.